

Comprehensive Characterization and Quantitation of Lipid Nanoparticle-Based RNA Formulations



Arthur Porfetye, Michelle Chen, Daniel Some — Wyatt Technology Corporation

Xiujuan Jia, Yong Liu, Angela Wagner, Yuejie Zhao, Katelyn J. Smith, Andreas M. Abend, Justin Pennington, Merck & CO., Inc.

Abstract

The therapeutic potential of lipid nanoparticles (LNPs) as delivery vehicles has been demonstrated in recent years cumulating in the current emergency use of the mRNA based SARS-CoV-2 vaccines. In order to ensure the safety and efficacy of the LNP-RNA vaccine or therapeutics, various quality attributes of LNP-RNA products need to be measured throughout the product development cycle. In this poster, we will demonstrate the use of a dynamics light scattering (DLS) Plate Reader for fast screening and quality control of the LNP preparations, and multi-angle light scattering (MALS) combined with ultraviolet (UV) and refractive index (dRI) detectors following size exclusion chromatography (SEC) or field flow fractionation (FFF) separation for in-depth characterization. SEC or FFF provides sized based separation, and MALS, UV, and dRI detectors enables biophysical characterization. We will show the data from the online detectors for the quantitation of particle size distribution, particle concentration, molecular weights of RNA and lipids, and sized-based RNA payload distribution of the LNP-RNA samples.

LNP-RNA Physical Attributes and Assays

Attribute	Assay	SEC/FFF-MALS-UV-dRI
mRNA integrity	Gel, qPCR	✓
LNP size	DLS analysis	✓
LNP distribution	DLS analysis	✓
Physical stability	DLS analysis	✓
LNP number	NTA analysis	✓
LNP morphology	TEM, cryo-EM	✓ (R_g/R_h)
LNP charge	Zeta potential	Possibly with EAF4
Encapsulation efficiency	Fluorescence	✓, new
mRNA concentration	Fluorescence	✓, new
Lipid concentration	LC-MS analysis	✓, new

Wyatt Solutions for Quantifying LNP-RNA Attributes

DLS Plate Reader

Size distribution and particle concentration with built-in automation

Mobius

Automated zeta potential measurement with an autosampler

ultraDAWN


Real-time LNP size and concentration for process development and PAT

SEC-MD


✓ MD (online multiple-detector): DAWN (MALS-DLS), UV (260 nm), Optilab (dRI)
 ✓ Software packages are 21 CFR Part 11 compliant, IQOQ service available

FFF-MD


Size-Based Separation by both SEC and FFF: mRNA and LNP




Comprehensive Characterization and Payload Analysis



Online LNP Payload Analysis: Cross Verification



Step-by-Step ASTRA LNP Payload Analysis Procedure



Conclusions

SEC-MALS/FFF-MALS Multi-Detector Analysis

LNP-RNA ATTRIBUTES QUANTIFIED

TECHNIQUES	Size Distribution	Lipid Concentration	LC-MS, LC-CAD, LC-ELSD
	Particle Number	RNA Concentration	LC-UV, Fluorescence
	Stability	Encapsulation Efficiency	
	Morphology	Size based Payload	Multiple assays with collected fractions

- SEC or FFF separates LNP with high resolution.
- Online detectors (MALS, UV at 260 nm, and dRI) provide comprehensive characterization and multi-attribute quantitation.
- The new LNP Analysis method enables size-based nucleic acid payload.
- MD-SEC and MD-FFF are essential tools for measuring LNP size, concentration, payload, and product quality. Software packages are 21 CFR 11 compliant, making both systems suitable in QC environment.
- Both systems are automated, robust, easy to adopt, with minimum hands-on time, less prone to experimental errors.
- MD-SEC and MD-FFF are orthogonal and complementary to the other conventional assays.